

WHICH VETTEC PRODUCTS ARE RIGHT FOR YOU?

PAD & PACKING MATERIALS

ADHESIVES



Equi-Pak

Soft & Supportive instant pad material

Specifications

46121

Available in:
160cc Cartridge

Color: Clear

Set Time:
30 seconds

Overview

Protect, cushion & support

Great for founder and laminitic cases

Bonds to sole & frog

Apply in thin protective layer or toward ground level for full support

No pad needed

Equi-Pak Soft

Extra soft instant pad material for use under pads

Specifications

46118

Available in:
160cc Cartridge

Color: Clear

Set Time:
40 seconds

Overview

Softer pad material

Can be used on horses with thin soles and deep commissures

Great for therapeutic cases

Designed for use under pads

Equi-Pak CS

Soft & supportive instant pad material for thrush

Specifications

46122

Available in:
160cc Cartridge

Color: Green

Set Time:
1 minute

Overview

Effectively deal with most cases of thrush

Provides durable, shock absorbing support

Eliminates needs to pick out feet or apply daily medication

Use as a preventive measure before thrush even begins

Sole-Guard

Sole protection & support for the unshod foot

Specifications

46160

Available in:
160cc Cartridge

Color: Black

Set Time:
30 seconds

Overview

Fast, effective protection and support for the sound footed, unshod horse

Sets firm to protect, but retains flexibility to move naturally with the foot

Provides protection and support for the bottom of the hoof in one easy application

Sil-Pak

Fast setting silicone packing material for use under pads

Specifications

46190

Available in:
160cc Cartridge

Color: Purple

Set Time:
1-2 minutes

Overview

Silicone material for use under pad in order to cushion and protect foot

Liquid formula fills all void areas minimizing debris migration providing a clean sole and eliminating odors

Super Fast

Fast setting hoof adhesive for foal extensions & custom shoes

Specifications

46140

Available in: 160cc & 50cc Cartridges

Color: Clear

Initial Set Time:
30 seconds

Weight bearing:
2-3 minutes

Overview

Ideal for foal extensions and making custom shoes

Great for building heels, correcting heel angles, and medial/lateral extensions

Good for small hoof repairs

Adhere

Glue-on shoe & hoof repair adhesive

Specifications

46150 & 46155

Available in: 160cc & 50cc Cartridges

Color: Black & Beige

Initial Set Time:
1 minute

Weight bearing:
4-5 minutes

Full cure: 6 minutes

Overview

For small and large hoof repairs and for gluing on shoes

Will bond to steel and aluminum

Excellent for urethane glue-on shoes and boots

Flexes like hoof wall and can be nailed